

# Mouse Brca1/Breast cancer type 1 susceptibility protein homolog ELISA Kit

Catalog No: #EK19391

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse Brca1/Breast cancer type 1 susceptibility protein homolog ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Breast cancer type 1 susceptibility protein homolog
Target Name	Breast cancer type 1 susceptibility protein homolog
Other Names	RING-type E3 ubiquitin transferase BRCA1
Accession No.	P48754
Uniprot	P48754
GeneID	12189;

## Application Details

Assay Type: Sandwich Detect Range: 78-5000pg/mL Sensitivity: 39.6pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: Brca1

Note: This product is for in vitro research use only