

Rat Brca1/Breast cancer type 1 susceptibility protein homolog ELISA Kit



Catalog No: #EK19392

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rat Brca1/Breast cancer type 1 susceptibility protein homolog ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Breast cancer type 1 susceptibility protein homolog
Target Name	Breast cancer type 1 susceptibility protein homolog
Other Names	RING-type E3 ubiquitin transferase BRCA1
Accession No.	O54952
Uniprot	O54952
GeneID	497672;

Application Details

Assay Type: Sandwich Detect Range: 0.78-50ng/mL Sensitivity: 0.401ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Brca1

Note: This product is for in vitro research use only