

# Human CD6/T-cell differentiation antigen CD6 ELISA Kit



Catalog No: #EK19393

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Human CD6/T-cell differentiation antigen CD6 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human T-cell differentiation antigen CD6
Target Name	T-cell differentiation antigen CD6
Other Names	T12, TP120, CD6
Accession No.	P30203
Uniprot	P30203
GeneID	923;

## Application Details

Assay Type: Sandwich Detect Range: 0.156-10ng/mL Sensitivity: 0.08ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: CD6

Note: This product is for in vitro research use only