

# Mouse Cd6/T-cell differentiation antigen CD6 ELISA Kit

Catalog No: #EK19394

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse Cd6/T-cell differentiation antigen CD6 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse T-cell differentiation antigen CD6
Target Name	T-cell differentiation antigen CD6
Other Names	CD6
Accession No.	Q61003
Uniprot	Q61003
GeneID	12511;

## Application Details

Assay Type: Sandwich Detect Range: 0.78-50ng/ml Sensitivity: 0.21ng/ml Sample Type: Serum, plasma, tissue homogenates and other biological fluids. Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: Cd6

Note: This product is for in vitro research use only