

Mouse Cd8a/T-cell surface glycoprotein CD8 alpha chain ELISA Kit

Catalog No: #EK19395

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Cd8a/T-cell surface glycoprotein CD8 alpha chain ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse T-cell surface glycoprotein CD8 alpha chain
Target Name	T-cell surface glycoprotein CD8 alpha chain
Other Names	T-cell surface glycoprotein Lyt-2,CD8a,Lyt-2,Lyt2
Accession No.	P01731
Uniprot	P01731
GeneID	12525;

Application Details

Assay Type: Sandwich Detect Range:0.312-20ng/ml Sensitivity:0.165ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name: Cd8a

Note: This product is for in vitro research use only