

Rat Cd8a/T-cell surface glycoprotein CD8 alpha chain ELISA Kit



Catalog No: #EK19396

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rat Cd8a/T-cell surface glycoprotein CD8 alpha chain ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat T-cell surface glycoprotein CD8 alpha chain
Target Name	T-cell surface glycoprotein CD8 alpha chain
Other Names	CD8 antigen 32 kDa chain,OX-8 membrane antigen,CD8a
Accession No.	P07725
Uniprot	P07725
GeneID	24930;

Application Details

Assay Type: Sandwich Detect Range:0.312-20ng/ml Sensitivity:0.101ng/ml Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name: Cd8a

Note: This product is for in vitro research use only