

Bovine ABCG2/ATP-binding cassette sub-family G member 2 ELISA Kit

Catalog No: #EK19397

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Bovine ABCG2/ATP-binding cassette sub-family G member 2 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Bovine
Specificity	Natural and recombinant Bovine ATP-binding cassette sub-family G member 2
Target Name	ATP-binding cassette sub-family G member 2
Other Names	Urate exporter, CD338
Accession No.	Q4GZT4
Uniprot	Q4GZT4
GeneID	536203;

Application Details

Assay Type: Sandwich Detect Range: 0.156-10ng/ml Sensitivity: 0.081ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: ABCG2

Note: This product is for in vitro research use only