

Mouse Abcg2/ATP-binding cassette sub-family G member 2 ELISA Kit



Catalog No: #EK19398

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Abcg2/ATP-binding cassette sub-family G member 2 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse ATP-binding cassette sub-family G member 2
Target Name	ATP-binding cassette sub-family G member 2
Other Names	Breast cancer resistance protein 1 homolog,Urate exporter,CD338,Abcp,Bcrp1
Accession No.	Q7TMS5
Uniprot	Q7TMS5
GeneID	26357;

Application Details

Assay Type: Sandwich Detect Range:0.156-10ng/mL Sensitivity:0.082ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Abcg2

Note: This product is for in vitro research use only