

# Pig ABCG2/ATP-binding cassette sub-family G member 2 ELISA Kit

Catalog No: #EK19399

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Pig ABCG2/ATP-binding cassette sub-family G member 2 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Pig
Specificity	Natural and recombinant Pig ATP-binding cassette sub-family G member 2
Target Name	ATP-binding cassette sub-family G member 2
Other Names	Brain multidrug resistance protein,Urate exporter,CD338,BMDP
Accession No.	Q8MIB3
Uniprot	Q8MIB3
GeneID	397073;

## Application Details

Assay Type: Sandwich Detect Range:1.56-100pg/mL Sensitivity:0.78pg/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene\_full\_name:ABCG2

Note: This product is for in vitro research use only