

Rat Abcg2/ATP-binding cassette sub-family G member 2 ELISA Kit

Catalog No: #EK19400

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rat Abcg2/ATP-binding cassette sub-family G member 2 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat ATP-binding cassette sub-family G member 2
Target Name	ATP-binding cassette sub-family G member 2
Other Names	Breast cancer resistance protein 1 homolog,Urate exporter,CD338,Bcrp1
Accession No.	Q80W57
Uniprot	Q80W57
GeneID	312382;

Application Details

Assay Type: Sandwich Detect Range:15.6-1000pg/mL Sensitivity:8.06pg/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Abcg2

Note: This product is for in vitro research use only