

## Rat Ambp/Protein AMBP ELISA Kit

Catalog No: #EK19401

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |  |
|--------------------|--|
| Product Name       | Rat Ambp/Protein AMBP ELISA Kit          |
| Brief Description  | ELISA Kits                               |
| Applications       | ELISA                                    |
| Species Reactivity | Rt                                       |
| Specificity        | Natural and recombinant Rat Protein AMBP |
| Target Name        | Protein AMBP                             |
| Other Names        | Itil                                     |
| Accession No.      | Q64240                                   |
| Uniprot            | Q64240                                   |
| GeneID             | 25377;                                   |

## Application Details

Assay Type: Sandwich Detect Range: 0.78-50ng/mL Sensitivity: 0.404ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: Ambp

Note: This product is for in vitro research use only