

Human AGTR2/Type-2 angiotensin II receptor ELISA Kit

Catalog No: #EK19402



Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Human AGTR2/Type-2 angiotensin II receptor ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Type-2 angiotensin II receptor
Target Name	Type-2 angiotensin II receptor
Other Names	Angiotensin II type-2 receptor,AT2
Accession No.	P50052
Uniprot	P50052
GeneID	186;

Application Details

Assay Type: SandwichDetect Range:1.56-100ng/mL Sensitivity:0.78ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nmgene_full_name:AGTR2

Note: This product is for in vitro research use only