

## Mouse Agtr2/Type-2 angiotensin II receptor ELISA Kit

Catalog No: #EK19403

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse Agtr2/Type-2 angiotensin II receptor ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Type-2 angiotensin II receptor
Target Name	Type-2 angiotensin II receptor
Other Names	Angiotensin II type-2 receptor,AT2
Accession No.	P35374
Uniprot	P35374
GeneID	11609;

## Application Details

Assay Type: Sandwich Detect Range:0.625-40ng/mL Sensitivity:0.292ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene\_full\_name:Agtr2

Note: This product is for in vitro research use only