

Mouse Cfi/Complement factor I ELISA Kit

Catalog No: #EK19404

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Cfi/Complement factor I ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Complement factor I
Target Name	Complement factor I
Other Names	C3B/C4B inactivator, If
Accession No.	Q61129
Uniprot	Q61129
GeneID	12630;

Application Details

Assay Type: Sandwich Detect Range: 0.312-20ng/mL Sensitivity: 0.158ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Cfi

Note: This product is for in vitro research use only