

## Rat Cfi/Complement factor I ELISA Kit

Catalog No: #EK19405

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Rat Cfi/Complement factor I ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Complement factor I
Target Name	Complement factor I
Other Names	C3B/C4B inactivator, If
Accession No.	Q9WUW3
Uniprot	Q9WUW3
GeneID	79126;

## Application Details

Assay Type: Sandwich Detect Range: 0.156-10ng/mL Sensitivity: 0.045ng/ml Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: Cfi

Note: This product is for in vitro research use only