

# Mouse Amacr/Alpha-methylacyl-CoA racemase ELISA Kit



Catalog No: #EK19406

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse Amacr/Alpha-methylacyl-CoA racemase ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Alpha-methylacyl-CoA racemase
Target Name	Alpha-methylacyl-CoA racemase
Other Names	2-methylacyl-CoA racemase,Macr1
Accession No.	O09174
Uniprot	O09174
GeneID	17117;

## Application Details

Assay Type: Sandwich Detect Range: 78-5000pg/mL Sensitivity: 39.2pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: Amacr

Note: This product is for in vitro research use only