

## Human ACP2/Lysosomal acid phosphatase ELISA Kit

Catalog No: #EK19407

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Human ACP2/Lysosomal acid phosphatase ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Lysosomal acid phosphatase
Target Name	Lysosomal acid phosphatase
Other Names	LAP
Accession No.	P11117
Uniprot	P11117
GeneID	53;

## Application Details

Assay Type: Sandwich Detect Range:0.156-10ng/mL Sensitivity:0.066ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene\_full\_name:ACP2

Note: This product is for in vitro research use only