

Human CD3G/T-cell surface glycoprotein CD3 gamma chain ELISA Kit

Catalog No: #EK19408

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Human CD3G/T-cell surface glycoprotein CD3 gamma chain ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human T-cell surface glycoprotein CD3 gamma chain
Target Name	T-cell surface glycoprotein CD3 gamma chain
Other Names	T-cell receptor T3 gamma chain, CD3g, T3G
Accession No.	P09693
Uniprot	P09693
GeneID	917;

Application Details

Assay Type: Sandwich Detect Range: 0.312-20ng/mL Sensitivity: 0.167ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: CD3G

Note: This product is for in vitro research use only