

Mouse Cd3g/T-cell surface glycoprotein CD3 gamma chain ELISA Kit

Catalog No: #EK19409

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Cd3g/T-cell surface glycoprotein CD3 gamma chain ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse T-cell surface glycoprotein CD3 gamma chain
Target Name	T-cell surface glycoprotein CD3 gamma chain
Other Names	T-cell receptor T3 gamma chain,CD3g
Accession No.	P11942
Uniprot	P11942
GeneID	12502;

Application Details

Assay Type: Sandwich Detect Range:0.156-10ng/mL Sensitivity:0.056ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name: Cd3g

Note: This product is for in vitro research use only