

Rat Cox4i1/Cytochrome c oxidase subunit 4 isoform 1, mitochondrial ELISA Kit



Catalog No: #EK19410

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rat Cox4i1/Cytochrome c oxidase subunit 4 isoform 1, mitochondrial ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Cytochrome c oxidase subunit 4 isoform 1, mitochondrial
Target Name	Cytochrome c oxidase subunit 4 isoform 1, mitochondrial
Other Names	Cytochrome c oxidase polypeptide IV, Cytochrome c oxidase subunit IV isoform 1, COX IV-1, Cox4, Cox4a
Accession No.	P10888
Uniprot	P10888
GeneID	29445;

Application Details

Assay Type: Sandwich Detect Range: 0.312-20ng/mL Sensitivity: 0.1ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Cox4i1

Note: This product is for in vitro research use only