

Mouse Cox5a/Cytochrome c oxidase subunit 5A, mitochondrial ELISA Kit



Catalog No: #EK19411

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Cox5a/Cytochrome c oxidase subunit 5A, mitochondrial ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Cytochrome c oxidase subunit 5A, mitochondrial
Target Name	Cytochrome c oxidase subunit 5A, mitochondrial
Other Names	Cytochrome c oxidase polypeptide Va
Accession No.	P12787
Uniprot	P12787
GeneID	12858;

Application Details

Assay Type: Sandwich Detect Range: 0.312-20ng/mL Sensitivity: 0.071ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Cox5a

Note: This product is for in vitro research use only