

OR13A1 Antibody

Catalog No: #44649

Package Size: #44649-1 50ul #44649-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	OR13A1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu
Specificity	OR13A1 Antibody detects endogenous levels of total OR13A1
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human OR13A1
Target Name	OR13A1
Other Names	Olfactory receptor 13A1
Accession No.	Swiss-Prot#: Q8NGR1/Q2M3M4NCBI Gene ID:79290/79290
Uniprot	Q8NGR1
GeneID	79290;
Calculated MW	36 kDa
Concentration	1.0mg/mL
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Images



Note: This product is for in vitro research use only