

Human CDH3/Cadherin-3 ELISA Kit

Catalog No: #EK19413

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Human CDH3/Cadherin-3 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Cadherin-3
Target Name	Cadherin-3
Other Names	Placental cadherin,P-cadherin,CDHP
Accession No.	P22223
Uniprot	P22223
GeneID	1001;

Application Details

Assay Type: Sandwich Detect Range:78-5000pg/mL Sensitivity:34pg/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids. Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:CDH3

Note: This product is for in vitro research use only