

Mouse Cd4/T-cell surface glycoprotein CD4 ELISA Kit

Catalog No: #EK19414

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Cd4/T-cell surface glycoprotein CD4 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse T-cell surface glycoprotein CD4
Target Name	T-cell surface glycoprotein CD4
Other Names	T-cell differentiation antigen L3T4, T-cell surface antigen T4/Leu-3, CD4
Accession No.	P06332
Uniprot	P06332
GeneID	12504;

Application Details

Assay Type: Sandwich Detect Range: 7.8-500ng/mL Sensitivity: 3.9ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids. Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Cd4

Note: This product is for in vitro research use only