

Rat Csnk2a1/Casein kinase II subunit alpha ELISA Kit

Catalog No: #EK19416

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|------------------------------------------------------------|
| Product Name | Rat Csnk2a1/Casein kinase II subunit alpha ELISA Kit |
| Brief Description | ELISA Kits |
| Applications | ELISA |
| Species Reactivity | Rt |
| Specificity | Natural and recombinant Rat Casein kinase II subunit alpha |
| Target Name | Casein kinase II subunit alpha |
| Other Names | CK II alpha |
| Accession No. | P19139 |
| Uniprot | P19139 |
| GeneID | 116549; |

Application Details

Assay Type: Sandwich Detect Range: 0.78-50ng/mL Sensitivity: 0.21ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Csnk2a1

Note: This product is for in vitro research use only