

## Human A2M/Alpha-2-macroglobulin ELISA Kit

Catalog No: #EK19417

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Human A2M/Alpha-2-macroglobulin ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Alpha-2-macroglobulin
Target Name	Alpha-2-macroglobulin
Other Names	Alpha-2-M,C3 and PZP-like alpha-2-macroglobulin domain-containing protein 5,CPAMD5,FWP007
Accession No.	P01023
Uniprot	P01023
GeneID	2;

## Application Details

Assay Type: Sandwich Detect Range:0.156-10ng/mL Sensitivity:0.078ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene\_full\_name:A2M

Note: This product is for in vitro research use only