

Mouse A2m/Alpha-2-macroglobulin ELISA Kit

Catalog No: #EK19418

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse A2m/Alpha-2-macroglobulin ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Alpha-2-macroglobulin
Target Name	Alpha-2-macroglobulin
Other Names	Alpha-2-M,Pregnancy zone protein,Pzp
Accession No.	Q61838
Uniprot	Q61838

Application Details

Assay Type: Sandwich Detect Range: 78-5000pg/ml Sensitivity: 39pg/ml Sample Type: Serum, plasma, tissue homogenates and other biological fluids. Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: A2m

Note: This product is for in vitro research use only