

Pig CBR1/Carbonyl reductase [NADPH] 1 ELISA Kit

Catalog No: #EK19419

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Pig CBR1/Carbonyl reductase [NADPH] 1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Pig
Specificity	Natural and recombinant Pig Carbonyl reductase [NADPH] 1
Target Name	Carbonyl reductase [NADPH] 1
Other Names	15-hydroxyprostaglandin dehydrogenase [NADP(+)],20-beta-hydroxysteroid dehydrogenase,NADPH-dependent carbonyl reductase 1,Prostaglandin 9-ketoreductase,Prostaglandin-E(2) 9-reductase,CBR,CRN
Accession No.	Q28960
Uniprot	Q28960
GeneID	397143;

Application Details

Assay Type: Sandwich Detect Range:0.156-10ng/mL Sensitivity:0.1ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:CBR1

Note: This product is for in vitro research use only