

## Mouse Ctnnb1/Catenin beta-1 ELISA Kit

Catalog No: #EK19420

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse Ctnnb1/Catenin beta-1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Catenin beta-1
Target Name	Catenin beta-1
Other Names	Beta-catenin, Catnb
Accession No.	Q02248
Uniprot	Q02248
GeneID	12387;

## Application Details

Assay Type: Sandwich Detect Range: 15.6-1000pg/mL Sensitivity: 7.85pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: Ctnnb1

Note: This product is for in vitro research use only