

Mouse Awms3/Uncharacterized protein C13orf33 homolog ELISA Kit

Catalog No: #EK19421

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Awms3/Uncharacterized protein C13orf33 homolog ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Uncharacterized protein C13orf33 homolog
Target Name	Uncharacterized protein C13orf33 homolog
Other Names	Activated in W/Wv mouse stomach 3,mAWMS3
Accession No.	Q14BA6
Uniprot	Q14BA6
GeneID	70717;

Application Details

Assay Type: Sandwich Detect Range: 78-5000pg/mL Sensitivity: 29.6 pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Awms3

Note: This product is for in vitro research use only