

Human BECN1/Beclin-1 ELISA Kit

Catalog No: #EK19422

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Human BECN1/Beclin-1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Beclin-1
Target Name	Beclin-1
Other Names	Coiled-coil myosin-like BCL2-interacting protein,Protein GT197,GT197
Accession No.	Q14457
Uniprot	Q14457
GeneID	8678;

Application Details

Assay Type: Sandwich Detect Range:78-5000pg/mL Sensitivity:38pg/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids. Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name: BECN1

Note: This product is for in vitro research use only