

## Chicken ALB/Serum albumin ELISA Kit

Catalog No: #EK19423

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Chicken ALB/Serum albumin ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Chicken
Specificity	Natural and recombinant Chicken Serum albumin
Target Name	Serum albumin
Other Names	Alpha-livetin
Accession No.	P19121
Uniprot	P19121

## Application Details

Assay Type: Sandwich
Detect Range:15.6-1000 pg/mL
Sensitivity:less than 7.85 pg/mL
Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul
Assay Time:3-5h
Detection wavelength:450 nm
gene_full_name:ALB

Note: This product is for in vitro research use only