

# Human ABCC3/Canalicular multispecific organic anion transporter 2 ELISA Kit



Catalog No: #EK19426

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |   |
|--------------------|---|
| Product Name       | Human ABCC3/Canalicular multispecific organic anion transporter 2 ELISA Kit   |
| Brief Description  | ELISA Kits  |
| Applications       | ELISA   |
| Species Reactivity | Hu  |
| Specificity        | Natural and recombinant Human Canalicular multispecific organic anion transporter 2   |
| Target Name        | Canalicular multispecific organic anion transporter 2   |
| Other Names        | ATP-binding cassette sub-family C member 3, Multi-specific organic anion transporter D, MOAT-D, Multidrug resistance-associated protein 3, CMOAT2, MLP2, MRP3 |
| Accession No.      | O15438  |
| Uniprot            | O15438  |
| GeneID             | 8714;   |

## Application Details

Assay Type: Sandwich Detect Range: 78-5000pg/mL Sensitivity: 33pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids. Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: ABCC3

Note: This product is for in vitro research use only