

# Rat Abcc6/Multidrug resistance-associated protein 6 ELISA Kit



Catalog No: #EK19429

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Rat Abcc6/Multidrug resistance-associated protein 6 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Multidrug resistance-associated protein 6
Target Name	Multidrug resistance-associated protein 6
Other Names	ATP-binding cassette sub-family C member 6,MRP-like protein 1,MLP-1,Mlp1,Mrp6
Accession No.	O88269
Uniprot	O88269
GeneID	81642;

## Application Details

Assay Type: Sandwich Detect Range:0.156-10ng/mL Sensitivity:0.088ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene\_full\_name:Abcc6

Note: This product is for in vitro research use only