

Mouse Abcc12/Multidrug resistance-associated protein 9 ELISA Kit

Catalog No: #EK19430

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Abcc12/Multidrug resistance-associated protein 9 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Multidrug resistance-associated protein 9
Target Name	Multidrug resistance-associated protein 9
Other Names	ATP-binding cassette sub-family C member 12,Mrp9
Accession No.	Q80WJ6
Uniprot	Q80WJ6
GeneID	244562;

Application Details

Assay Type: SandwichDetect Range:0.156-10ng/mL Sensitivity:0.061ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nmgene_full_name:Abcc12

Note: This product is for in vitro research use only