

# Rat Abcc12/Multidrug resistance-associated protein 9 ELISA Kit



Catalog No: #EK19431

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Rat Abcc12/Multidrug resistance-associated protein 9 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Multidrug resistance-associated protein 9
Target Name	Multidrug resistance-associated protein 9
Other Names	ATP-binding cassette sub-family C member 12,Mrp9
Accession No.	Q6Y306
Uniprot	Q6Y306
GeneID	291923;

## Application Details

Assay Type: SandwichDetect Range:78-5000pg/mL Sensitivity:32pg/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids. Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nmgene\_full\_name:Abcc12

Note: This product is for in vitro research use only