

Mouse Abcb4/Multidrug resistance protein 2 ELISA Kit

Catalog No: #EK19432

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Abcb4/Multidrug resistance protein 2 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Multidrug resistance protein 2
Target Name	Multidrug resistance protein 2
Other Names	ATP-binding cassette sub-family B member 4,P-glycoprotein 2,Mdr2,Pgy-2,Pgy2
Accession No.	P21440
Uniprot	P21440
GeneID	18670;

Application Details

Assay Type: Sandwich Detect Range:78-5000pg/mL Sensitivity:40.1pg/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Abcb4

Note: This product is for in vitro research use only