

OR2L2 Antibody

Catalog No: #44651

Package Size: #44651-1 50ul #44651-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	OR2L2 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu
Specificity	OR2L2 Antibody detects endogenous levels of total OR2L2
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human OR2L2
Target Name	OR2L2
Other Names	OR2L2, HSHTPCRH07, HTPCRH07, OR2L12, Olfactory receptor 2L12, Olfactory receptor 2L2, OR2L4P, Olfactory receptor 2L4, OR1-48
Accession No.	Swiss-Prot#: Q8NH16NCBI Gene ID:26246
Uniprot	Q8NH16
GeneID	26246;
Calculated MW	35 kDa
Concentration	1.0mg/mL
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Images



Note: This product is for in vitro research use only