

# Rat Abcb4/Phosphatidylcholine translocator ABCB4 ELISA Kit



Catalog No: #EK19433

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Rat Abcb4/Phosphatidylcholine translocator ABCB4 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Phosphatidylcholine translocator ABCB4
Target Name	Phosphatidylcholine translocator ABCB4
Other Names	ATP-binding cassette sub-family B member 4, Multidrug resistance protein 2, Multidrug resistance protein 3, P-glycoprotein 2, P-glycoprotein 3, Mdr2, Pgp3, Pgy2
Accession No.	Q08201
Uniprot	Q08201
GeneID	24891;

## Application Details

Assay Type: Sandwich Detect Range: 0.312-20ng/mL Sensitivity: 0.156ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: Abcb4

Note: This product is for in vitro research use only