

Rat Abcb11/Bile salt export pump ELISA Kit

Catalog No: #EK19434

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rat Abcb11/Bile salt export pump ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Bile salt export pump
Target Name	Bile salt export pump
Other Names	ATP-binding cassette sub-family B member 11,Sister of P-glycoprotein,Bsep,Spgp
Accession No.	O70127
Uniprot	O70127
GeneID	83569;

Application Details

Assay Type: Sandwich Detect Range:0.312-20ng/mL Sensitivity:0.07ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Abcb11

Note: This product is for in vitro research use only