

Human ACADM/Medium-chain specific acyl-CoA dehydrogenase, mitochondrial ELISA Kit



Catalog No: #EK19435

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Human ACADM/Medium-chain specific acyl-CoA dehydrogenase, mitochondrial ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Medium-chain specific acyl-CoA dehydrogenase, mitochondrial
Target Name	Medium-chain specific acyl-CoA dehydrogenase, mitochondrial
Other Names	MCAD
Accession No.	P11310
Uniprot	P11310
GeneID	34;

Application Details

Assay Type: Sandwich Detect Range: 0.156-10ng/mL Sensitivity: 0.078ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: ACADM

Note: This product is for in vitro research use only