

Mouse Acadm/Medium-chain specific acyl-CoA dehydrogenase, mitochondrial ELISA Kit

Catalog No: #EK19436

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Acadm/Medium-chain specific acyl-CoA dehydrogenase, mitochondrial ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Medium-chain specific acyl-CoA dehydrogenase, mitochondrial
Target Name	Medium-chain specific acyl-CoA dehydrogenase, mitochondrial
Other Names	MCAD
Accession No.	P45952
Uniprot	P45952
GeneID	11364;

Application Details

Assay Type: Sandwich Detect Range: 31.2-2000pg/mL Sensitivity: 15pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Acadm

Note: This product is for in vitro research use only