



Rat Cckbr/Gastrin/cholecystokinin type B receptor
ELISA Kit

Catalog No: #EK19439

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Rat Cckbr/Gastrin/cholecystokinin type B receptor ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Gastrin/cholecystokinin type B receptor
Target Name	Gastrin/cholecystokinin type B receptor
Other Names	CCK-B receptor,Cholecystokinin-2 receptor,CCK2-R
Accession No.	P30553
Uniprot	P30553
GeneID	25706;

Application Details

Assay Type: SandwichDetect Range:78-5000pg/ml Sensitivity:40.4pg/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nmgene_full_name:Cckbr

Note: This product is for in vitro research use only