

Mouse Apaf1/Apoptotic protease-activating factor 1 ELISA Kit

Catalog No: #EK19440

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Apaf1/Apoptotic protease-activating factor 1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Apoptotic protease-activating factor 1
Target Name	Apoptotic protease-activating factor 1
Other Names	APAF-1
Accession No.	O88879
Uniprot	O88879
GeneID	11783;

Application Details

Assay Type: Sandwich Detect Range:0.156-10ng/mL Sensitivity:0.078ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Apaf1

Note: This product is for in vitro research use only