

Human AIFM1/Apoptosis-inducing factor 1, mitochondrial ELISA Kit



Catalog No: #EK19441

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Human AIFM1/Apoptosis-inducing factor 1, mitochondrial ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Apoptosis-inducing factor 1, mitochondrial
Target Name	Apoptosis-inducing factor 1, mitochondrial
Other Names	Programmed cell death protein 8,AIF,PDCD8
Accession No.	O95831
Uniprot	O95831
GeneID	9131;

Application Details

Assay Type: Sandwich Detect Range:0.156-10ng/mL Sensitivity:0.06ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:AIFM1

Note: This product is for in vitro research use only