

Mouse Aifm1/Apoptosis-inducing factor 1, mitochondrial ELISA Kit



Catalog No: #EK19442

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Aifm1/Apoptosis-inducing factor 1, mitochondrial ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Apoptosis-inducing factor 1, mitochondrial
Target Name	Apoptosis-inducing factor 1, mitochondrial
Other Names	Programmed cell death protein 8,Aif,Pdcd8
Accession No.	Q9Z0X1
Uniprot	Q9Z0X1
GeneID	26926;

Application Details

Assay Type: Sandwich Detect Range:78-5000pg/mL Sensitivity:32pg/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids. Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Aifm1

Note: This product is for in vitro research use only