

Mouse Adrb1/Beta-1 adrenergic receptor ELISA Kit

Catalog No: #EK19444

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|--|
| Product Name | Mouse Adrb1/Beta-1 adrenergic receptor ELISA Kit |
| Brief Description | ELISA Kits |
| Applications | ELISA |
| Species Reactivity | Ms |
| Specificity | Natural and recombinant Mouse Beta-1 adrenergic receptor |
| Target Name | Beta-1 adrenergic receptor |
| Other Names | Beta-1 adrenoreceptor, Beta-1 adrenoceptor, Adrb1r |
| Accession No. | P34971 |
| Uniprot | P34971 |
| GeneID | 11554; |

Application Details

Assay Type: Sandwich Detect Range: 0.78-50ng/mL Sensitivity: 0.17ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Adrb1

Note: This product is for in vitro research use only