

Rat Cnr2/Cannabinoid receptor 2 ELISA Kit

Catalog No: #EK19446

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rat Cnr2/Cannabinoid receptor 2 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Cannabinoid receptor 2
Target Name	Cannabinoid receptor 2
Other Names	CB-2
Accession No.	Q9QZN9
Uniprot	Q9QZN9
GeneID	57302;

Application Details

Assay Type: Sandwich Detect Range: 1.56-100ng/mL Sensitivity: 0.83ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Cnr2

Note: This product is for in vitro research use only