

Bovine A1BG/Alpha-1B-glycoprotein ELISA Kit

Catalog No: #EK19447

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Bovine A1BG/Alpha-1B-glycoprotein ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Bovine
Specificity	Natural and recombinant Bovine Alpha-1B-glycoprotein
Target Name	Alpha-1B-glycoprotein
Other Names	Alpha-1-B glycoprotein
Accession No.	Q2KJF1
Uniprot	Q2KJF1
GeneID	518955;

Application Details

Assay Type: Sandwich Detect Range: 3.12-200ng/mL Sensitivity: 1.75ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: A1BG

Note: This product is for in vitro research use only