

Rat A1bg/Alpha-1B-glycoprotein ELISA Kit

Catalog No: #EK19450

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rat A1bg/Alpha-1B-glycoprotein ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Alpha-1B-glycoprotein
Target Name	Alpha-1B-glycoprotein
Other Names	Alpha-1-B glycoprotein,C44,Liver regeneration-related protein 1
Accession No.	Q9EPH1
Uniprot	Q9EPH1
GeneID	140656;

Application Details

Assay Type: Sandwich Detect Range:78-5000pg/mL Sensitivity:36pg/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids. Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:A1bg

Note: This product is for in vitro research use only